


STRAW-b bootcamp

Overview



STRAW, STRAW-b, P-ONE

km³-scale Detectors for Astrophysical Neutrinos

Point-source search with astrophysical neutrinos



The directions with highest sensitivity for astrophysical neutrinos per detector



P-ONE

- 2660m deep
- Use of the existing deep sea infrastructure of ONC

Our Path to a km³-scale Neutrino Detector

2018

STRAW

2020

STRAW-b





202X

P-ONE 10 String

202X

P-ONE

STRings for the Absorption-length in Water




Motivation & Design


STRAWb

STRAW-b - Sensors


Muon-Tracker




2x LiDAR




Standard + Art




WOM




Mini-Spectrometer





2x PMT-Spectrometer



Standard + Art




Standard + Art




STRAW-b - Sensors

 Minispectrometer
  Camera
  LED (Lucifer)


Muon-Tracker




2x LiDAR




Standard + Art




WOM




Mini-Spectrometer





2x PMT-Spectrometer



Standard + Art




Standard + Art



Motivation for STRAWb

 Minispectrometer
  Camera
  LED (Lucifer)

Module	Goal
Muon-Tracker	Measure Muons - reconstruct the direction & energy
PMT-Spectrometer	Bioluminescence/animal classification + spectra + image
PMT-Spectrometer + LiDAR	Wavelength dependent attenuation length of water
PMT-Spectrometer	Measure background characteristics for PONE (in addition to STRAW)
LiDAR	Attenuation length and scattering length of water
Mini/PMT-Spectrometer	Monitor sedimentation of module with camera
Logger	Analyse stress on construction during deployment
...	...




Status
STRAWb


Status

Minispectrometer Camera LED (Lucifer)


Muon-Tracker




2x LiDAR




Standard + Art




WOM




Mini-Spectrometer





2x PMT-Spectrometer



Standard + Art



Standard + Art




Status



▼ Minispectrometer
 ▼ Camera
 ▼ LED (Lucifer)

No connection since 28.11.2021
Investigation ongoing


Muon-Tracker




2x LiDAR




Mini-Spectrometer




2x PMT-Spectrometer



Standard + Art



Standard + Art



DAQ
STRAWb


DAQ

SDAQ - "Simple Data Acquisition"

Standalone python package:

<https://gitlab.lrz.de/ecp/strawb/mctl/-/tree/master/h5daq>


Sensor A



SDAQ Job

- Buffer for settings and data

Sensor B



SDAQ Job

- Buffer for settings and data

SDAQ Daemon

- Collects Data
- Does rollovers
- Compression

HDF5 File

- A hdf5-group per job

Group A

- Dataset Time [1D array]
- X Dataset(s) [nD array]
- First dimension is the time

Group B

- Dataset Time [1D array]
- X Dataset(s) [nD array]
- First dimension is the time

File Types

Module

.hdf5

- Voltage+Current readings per sensors
- Accel. Sensor
- Magnetic Sensor
- Temperature Sensors
- LED settings

Data-Logger

.hdf5

- Accel. Sensor
- Magnetic Sensor
- Temperature Sensors

- Connection works
- Accessing data fails

Camera

.hdf5

- Raw RGB-Intensity data
- Exposure time
- Gain settings
- Power settings

Standard Module

Specialised Module

Mini-Spectrometer

.hdf5

- Raw Wavelength/Intensity data
- Exposure time
- Sensor temperature
- Power settings

LiDAR

.hdf5

- Gimbal parameters
- PMT settings
- Laser settings
- Power settings

.hld / .txt

- Raw TOT data (Time over Threshold)

PMT-Spectrometer

.hdf5

- PMT rates with 1 kHz readout
- HV settings
- Threshold settings
- Power settings
- Raw TOT data

.hld / .txt

- Raw TOT data (Time over Threshold)

MuonTracker

.hdf5

- SiPM rates with 10 Hz readout
- HV settings
- Threshold settings
- Power settings
- Raw TOT data


.hld / .txt

- Raw TOT data (Time over Threshold)


Digital infrastructure

STRAWb


Digital Infrastructure for STRAWb



Module Monitoring



LRZ Setup



LRZ DSS Structure

DSS Container


- Mounted at the VM under: /dss/
- STRAW and STRAWb data



STRAWb Data

- All raw module data are located at: /dss/strawb/raw_module_data/
- Files are sorted with: <device_code>/<year-month>/
I.e.: tumlidar001/2021_10/...
- Includes a DB file:
../raw_module_data/pandas_file_sync_db.gz

A simple and reliable way to get the files.
No need to deal with path and filenames.



Example Notebook to generate the plot:

https://github.com/pone-software/strawb_package/blob/master/examples/pandas_file_sync_db_explorer.ipynb

Working with the Data

Python STRAWb package

Goal's

- Simplifies to access the data from the ONC server
- Simplifies reading of all STRAWb file types
- Collection of tools to process the data, including calibration
- Includes various examples for an easy start
- Features a DB for simple file access

Code

https://github.com/pone-software/strawb_package

The screenshot shows the GitHub repository page for 'pone-software/strawb_package'. The repository is private and has 3 watchers. The main branch is 'master', with 2 other branches and 0 tags. The repository has 209 commits. The commit history is as follows:





Commit	Message	Time
f04e9d8	Merge branch 'minispec_2nd_update' into camera_dev	last month
	added LiDAR to todo list	24 days ago
	fixed bug	23 days ago
	added docs to onc_vm_downloader.py	29 days ago
	updated imports	23 days ago
	added LiDAR test	24 days ago
	second main chance, simplify the functions, add examples	3 months ago
	Ignore the output of jupyter notebooks	28 days ago
	added `**/__pycache__`	last month
	Multiple improvements	5 months ago
	typo	23 days ago
	fixed config	29 days ago
	Multiple improvements	5 months ago
	added plotly, statsmodels	28 days ago
	Multiple improvements	5 months ago

PMT-Spectrometer



Motivation - Goals



Module	Goal
Muon-Tracker	Measure Muons - reconstruct the direction & energy
PMT-Spectrometer	Bioluminescence/animal classification + spectra + image
PMT-Spectrometer + LiDAR	Wavelength dependent attenuation length of water
PMT-Spectrometer	Measure background characteristics for PONE (in addition to STRAW)
LiDAR	Attenuation length and scattering length of water
Mini/PMT-Spectrometer	Monitor sedimentation of module with camera
Logger	Analyse stress on construction during deployment



Wavelength dependent attenuation


Measurement principle

- LED on the side and bottom of the LiDAR emit light
- Light propagates 24 m through water
- Measure the spectra with the PMTs, Mini-Spectrometer and Camera
- + calibration lab measurements, e.g.:



Not 100% clear if it works as only one distance is available
But it can be use as an in-situ calibration of the PMT spectra

PMT-Spectrometer - Design




PMT-Spectrometer - Optics

Optics are optimised for

- all 12 PMTs have the same field of view
- have the same acceptance in the field of view
- correct the optical effect of the spherical glass

The 12 PMT Detection Unit

- A PMT Detection Unit is the PMT, lens and filter
- The PMT detection units are symmetrical placed around the vertical axis in 3 groups
- Within each group the Detection Units have the same radius and the same declination angle to the vertical axis



- Inner Group
- Middle Group
- Outer Group

